

Assay Name: Cell cycle analysis using DAPI and BrdU

Assay ID: Celigo_02_0001

Description: Analysis of cell cycle using DAPI and BrdU incorporation on A549 cells

Stains: Alexa Fluor® 488 goat anti-mouse antibody and DAPI

Imaging channels: Green and Blue

Image analysis algorithm: Celigo software Expression Analysis: Target 1 + Mask


Methods:

1. Seed cells in 96-well plate and incubate cells overnight
2. Add 100 μ l of compound at 2x desired final concentration (e.g. Aphidicolin, an inhibitor of DNA synthesis)
3. BrdU incorporation step, 100 μ M BrdU and incubate for 2 hours
4. Fixation and DNA denaturation step, 200 μ l/well of Ethanol/Glycine solution for 30 minutes
5. Antibody labeling steps:
 - a. 50 μ l Anti-BrdU mouse monoclonal antibody and incubate for 2 hr
 - b. 50 μ l Anti-BrdU Alexa Fluor® 488 goat anti-mouse antibody and incubate for 1 hr and wash
6. Add 100 μ l DAPI mix to each well and incubate for 30 min
7. Image and gate on the Celigo


Results:

Gating interface displaying the horseshoe plot and cells are color-coded dependent on the phase of the cell cycle they are in.


Gating Interface Plot



Cell Segmentation




- BrdU integrated intensities vs DAPI integrated intensities are plotted to create the horseshoe shape in the above left image.
- S-phase cells are identified in green, G0/G1 cells in purple and G2/M in blue in a control well in the above right image.



Aphidicolin Conc. Response

- 96 well plate readout displaying % cells in S-phase and a heat map display.
- The top half of the plate contains 100nM Aphidicolin and controls; the bottom is a serial dilution of Aphidicolin across the plate.
- In the heat map, green represents low % cells in S-Phase, yellow and intermediate % cells in S-phase and red is a high % of cells in S-phase.

Effect of Aphidicolin on Cell Cycle



	G0/G1	S-Phase	G2/M
IC50	2.11e-007	2.952e-007	7.449e-007