

## Assay Name: Transwell migration using DAPI

Assay ID: Celigo\_04\_0001



### Table of Contents

Experiment: Identification of cells which have migrated through the membrane of a transwell insert.....	2
Celigo Setup.....	2
Assay Protocol and Plate Setup.....	2
Results .....	3
1. Bright field and Blue image of well A1 .....	3
2. Numerical counts of DAPI positive .....	4
Conclusion .....	5

## Experiment: Transwell migration using DAPI

Purpose	Image and count the number of cells which have migrated through the membrane of a transwell insert using fluorescent detection of DAPI-stained cells.
Current Method(s)	Manual visualization using standard light microscope
Target Cell Type	Three different cancer cell lines (undisclosed)
Experiment Plan	Scan plate using Bright field and Blue channels
Hypothesis	We will be able to stain and count the number of cells which have migrated through the membrane of a transwell insert.

## Celigo Setup

Plate Type	Transwell inserts placed in a 24-well BD Falcon 353047 plate
Scan Channels	Bright field and Blue channels
Resolution	1 $\mu\text{m}$ /pixel
Scan Area	Partial Well - 25 FOV
Analysis Method	Target 1 + 2
Scan Frequency	End Point
Scan Time	3 minutes



## Assay Protocol and Plate Setup

1. Three different cancer cell lines were seeded into a transwell insert and grown according to the customer's standard procedure, with or without treatment to stimulate migration.
2. 24hrs after seeding, the transwell inserts were removed and the cells were fixed in the transwell insert with formaldehyde for 10 minutes.
3. Cells were then washed with water to remove the formaldehyde.
4. Using a sterile cotton swap, cells which had not migrated through the membrane were scraped off the top of the transwell insert.
5. A staining solution of DAPI, with a permeability reagent, Triton X-100, was prepared in 5 mL PBS.
  - a. 1:1,000 dilution of DAPI stock (10mg/mL) for a final 10  $\mu\text{g}/\text{mL}$ 
    - i. Nexcelom, Cat# CS1-0127
  - b. Add Triton X-100 for a final 1% solution
6. 600  $\mu\text{L}$  of staining solution was added to a clean well of a 24-well BD Falcon plate (Cat# 353047).
7. Placed transwell insert into a well with staining solution and allowed to stain for approximately 10-15 minutes.
8. Removed transwell insert and washed with PBS to remove any stain or Triton X-100.
9. Added 600  $\mu\text{L}$  of PBS to a clean well of a 24-well BD Falcon plate.
10. Placed washed transwell insert into the well.
11. Imaged on Celigo using the Target 1 + 2 application with Bright field and Blue channels.



## Results

1. Bright field and Blue image of well A1

Full insert images / partial well images (25 FOV)




Zoomed-in images




Note: The circular objects seen in the Bright field images are the pores in the transwell membrane. The red arrow points to one of the many pores that can be seen in the image.

## Overlay of Bright field and Blue images



**Note:** The overlay of the bright field image with the DAPI image allows the confirmation that the circular objects seen in the bright field image are the pores in the membrane, and not cells (see red arrow). The cells are stained with DAPI and appear as the blue dots in the image (see white arrow). Some cells are still migrating through the pores, thus the blue dots occasionally line up with some of the pores in the image (see yellow arrow).

## 2. Numerical counts of DAPI-positive cells



## Conclusion

- The Celigo successfully imaged and counted the number of DAPI-positive stained cells which had migrated through a membrane in a transwell insert.
- Acquisition of high resolution DAPI and bright field images of the transwell membrane surface in a 24 well plate format took ~ 3 minutes.