

1. ABSTRACT

Apoptosis and cell cycle play an important role in various aspects of preclinical pharmaceutical drug discovery and validation. The ability to quickly determine the cytotoxic effect of chemical compounds on cancer cells allows researchers to efficiently identify potential drug candidates for further development in the pharmaceutical discovery pipeline. Recently, a plate-based imaging cytometry system, Celigo Imaging Cytometer, has been used to for high-throughput fluorescence cell cycle and apoptosis analysis. In this study, we demonstrate the use of Celigo imaging cytometry for apoptosis and cell cycle detection by studying the dose response effect of nocodazole on cell cycle and staurosporine on apoptosis. For cell cycle analysis, the cells are labeled with propidium iodide and BrdU. For apoptosis analysis, the cells are labeled with Annexin V-PE and Hoechst 33342. The experimental results were evaluated to validate the imaging cytometric capabilities of the Imaging Cytometry system. The plate-based imaging cytometer utilizes bright-field and three fluorescence channels (Blue, Green, and Red) for multi-channel analysis. By utilizing the F theta lens technology, uniform bright-field image is captured for more accurate cell counting and analysis of the entire well. In addition, Celigo analysis software is used to report numerous parameters allowing detailed fluorescence-based cell population characterization. The ability of Celigo to rapidly and cost-effectively perform plate-based fluorescent assays has the potential of improving research efficiency, especially for adherent cells where plate-based cytometer does not require trypsinization for cell population analysis.



2. CELIGO IMAGING CYTOMETRY FOR APOPTOSIS AND CELL CYCLE ANALYSIS

1. Celigo Imaging Cytometer is a plate-based cytometer that can scan the entire well of standard microplates and captures bright-field and fluorescent images

2. The captured images are analyzed with the Celigo software to measure size, morphology, cell count, confluence, and fluorescent intensity

3. The measured parameters are used to generate cell proliferation kinetic data, GFP/RFP expression, tumor spheroid size change, DNA cell cycle analysis, apoptosis, and ADCC cytotoxicity results

3. CELL CYCLE ANALYSIS USING CELIGO IMAGE CYTOMETRY




- Cell cycle analysis is often used in drug screening to identify compounds that can affect cell proliferation and growth
- Celigo image cytometer can be used to measure the G0/G1, S, and G2 phase in a cell population to determine the effect of treatment
- In this work, A549 cells are stained with DAPI, a fluorescent nuclear dye, and BrdU, a thymidine analog, to incorporate into cells for detection of DNA content and replication
- BrdU uses a primary anti-mouse monoclonal antibody and anti-Gam AF-488 secondary antibody
- A549 cells are treated with nocodazole for 18 hours, at 150 nM (middle) and 50 μM (right)
 - The results showed reduction in G0/G1 phase at 150 nM, and increase in apoptotic cells at 50 μM
- A549 cells are treated with aphidicolin for 18 hours, from 1 nM (left) to 100 μM (middle)
 - The results showed large S phase population, while at high dosage, that population is completely eliminated
- The aphidicolin EC50 for G0/G1 and S phases, are 463 and 335 nM, respectively

6. SUMMARY AND CONCLUSION

In conclusion, the Celigo image cytometer has been demonstrated to perform cell viability, cell cycle, and apoptosis analysis, which are three critical cell-based assays for immunological, oncological, and toxicological research. By utilizing an automated image cytometry system, drug screen results can be quickly obtained using high throughput plate format, and significantly increase the productivity of research capabilities.

References
 • Sasaki et al. Flow cytometric estimation of cell cycle parameters using a monoclonal antibody to bromodeoxyuridine. *Cytometry*. 1986; 7(4):391-5.

4. CELL VIABILITY ANALYSIS USING CELIGO IMAGE CYTOMETRY



5. APOPTOSIS ANALYSIS USING CELIGO IMAGE CYTOMETRY

