



Novel image cytometric method for detection of physiological and metabolic changes in *Saccharomyces cerevisiae*

Leo L. Chan^{1,2}, Alexandria L. Kury^{1,2}, Alisha Wilkinson^{1,2}, Charlotte Berkes², Alnoor Pirani¹
¹Nexcelom Bioscience LLC, 360 Merrimack St. Building 9, Lawrence, MA 01843
²Center for Biotechnology and Biomedical Sciences, Merrimack College, North Andover, MA 01845

MBAA Annual Conference
 October 23-26, 2013
 Hilton Tucson Austin
 Austin, TX


1. ABSTRACT

The study and monitoring of physiological and metabolic changes in *Saccharomyces cerevisiae* (*S. cerevisiae*) has been a key research area for the brewing, baking and biofuels industries, which rely on these economically important yeasts to produce their products. Specifically for breweries, physiological and metabolic parameters such as viability, vitality, glycogen, neutral lipid, and trehalose content can be measured to better understand the status of *S. cerevisiae* during fermentation. Traditionally, these physiological and metabolic changes can be qualitatively observed using fluorescence microscopy or flow cytometry for quantitative fluorescence analysis of fluorescently labeled cellular components associated with each parameter. However, both methods pose known challenges to the end-users. Specifically, conventional fluorescent microscopes lack automation and fluorescence analysis capabilities to quantitatively analyze large numbers of cells. Although flow cytometry is suitable for quantitative analysis of tens of thousands of fluorescently labeled cells, the instruments require a considerable amount of maintenance, highly trained technicians, and the system is relatively expensive to both purchase and maintain. In this work, we demonstrate the first use of Cellometer Vision for the kinetic detection and analysis of vitality, glycogen, neutral lipid, and trehalose content of *S. cerevisiae*. This method provides an important research tool for large and small breweries to study and monitor these physiological behaviors during production, which can improve fermentation conditions to produce consistent and higher quality products.


2. CURRENT METHODS FOR MEASURING METABOLIC CHANGES IN YEAST

Methods	Description	Known Issues
Hemocytometer	• Manually counting viable cells	• Time-consuming and tedious process • Requires experienced user for accurate counting
Fluorescence Microscopy	• Visualization of fluorescently-stained yeast cells	• Qualitative observe instead of quantitative analysis • Not automated, low throughput
Flow Cytometry	• Quantitative fluorescent analysis • Automated analysis	• Relatively expensive and high maintenance • Requires experienced user for proper operation • Cannot visually observe budding yeasts

3. CELLOMETER IMAGE CYTOMETRY INSTRUMENTATION



4. FLUORESCENCE INTENSITY ANALYSIS USING IMAGE CYTOMETRY



- Target yeast samples are stained with the appropriate fluorescent labels
- Using Cellometer image cytometer to capture bright-field and fluorescent images
- The Cellometer software counts all the yeast cells and outlines the shape of the yeast in bright-field
- The mask that is generated is used in fluorescent images to measure intensity from each yeast particle
- The fluorescence intensity data is then exported to FCS Express 4 for data analysis and presentation


5. DESCRIPTION OF PHYSIOLOGICAL AND METABOLIC PARAMETERS

Parameters	Description
Viability	• Measurement of percent live cells in the sample • Determine amount of viable yeast used for fermentation • Improve consistency by utilizing consistent amount of live cells
Vitality	• Measurement of percent of cells able to perform enzymatic reactions • Metabolically active cells are fermenting during production • Vitality and viability do not usually correlate
Trehalose	• Sugar disaccharide that supplies energy during cell cycle • Protects against stress, high EtOH, heat, dehydration, oxidation, pH • Mobilized quickly, related to survival capacity
Neutral Lipid	• Energy-rich molecules that are stored in yeasts • Cannot be metabolized • Protect from toxic substance like EtOH, prolong survival capacity
Glycogen	• Macromolecule that provides the energy and carbohydrates required for yeast sterols and lipid synthesis • Mobilized quickly to provide energy during and after fermentation • It has been shown to correlate to vitality

6. DETECTION METHOD OF PHYSIOLOGICAL AND METABOLIC PARAMETERS


Parameters	Detection Method
Viability	• Measured by Acridine Orange (AO, EX: 475 nm, EM: 535 nm) and Propidium Iodide (PI, EX: 525 nm, EM: 595 nm) fluorescent nucleic acid staining method
Vitality	• Measured by Carboxyfluorescein Diacetate (CFDA) fluorescent enzymatic staining method (EX: 475 nm, EM: 535 nm)
Trehalose	• Measured by staining with Concanavalin A-FITC (EX: 475 nm, EM: 535 nm)
Neutral Lipid	• Measured by staining with Nile Red (EX: 525 nm, EM: 600 nm)
Glycogen	• Measured by staining with Acriflavine (EX: 475 nm, EM: 535 nm)

7. EXPERIMENTAL PROTOCOL FOR BUDDING MEASUREMENT



- Yeast colonies are streaked out on an agar plate
- A 5 ml YPD test tube is inoculated and allowed to culture overnight
- A stationary phase sample of 800 µl is transferred to a tube of 20 ml YPD media
 - Samples are collected at 5, 10, 25, and 30 hours
- Collected yeast samples are resuspended in PBS or EtOH
- For PBS: stained with AO/PI, CFDA; for EtOH: stained with Acriflavine, Nile Red, Con A-FITC

8. PHYSIOLOGICAL AND METABOLIC PARAMETER TIME-COURSE MEASUREMENT



9. TIME-COURSE MEASUREMENT RESULTS AND CONCLUSION

