

A Rapid Alternative Method for Cell Cycle Analysis Using Cellometer Vision

Leo L. Chan¹, Ph.D., Xuemei Zhong², Ph.D., Jean Qiu¹, Ph.D., Peter Y. Li¹, Ph.D., and Bo Lin¹, Ph.D.
¹Nexcelom Bioscience LLC, 360 Merrimack St. Building 9, Lawrence, MA 01843
²Hematology Oncology Section, Department of Medicine, Boston University Medical Center, Boston, MA 02118

1. ABSTRACT

Traditional cell cycle analysis instrumentation involves the use of fluorescence microscopy, laser scanning cytometry, and flow cytometry. In this work, we demonstrate new applications for the Cellometer Vision, which utilizes an imaging cytometry method for cell cycle analysis. This method employs both bright-field and fluorescence imaging of a disposable cell counting chamber to quickly provide concentration and percentages of cell subpopulations. Signal measurement and data analysis could be performed less than five min. Cell cycle assay was demonstrated by detection of changes in phase percentage within a cell population when Jurkat cells were induced with three different drugs, nocodazole, etoposide, and aphidicolin that arrest cell cycle at G₂/M, S, and G₀/G₁ phase, respectively. Jurkat cells were fixed with 75% ethanol and stained with propidium iodide. Individual cells were identified by the software and the fluorescence intensity was measured from within the cell. The imaging cytometry method was able to detect and analyze fluorescence from the Jurkat cells, which generated cell cycle data plot for each drug at three different concentrations. The cell cycle analysis results were highly consistent with the results obtained from FACS-calibur flow cytometer. This imaging cytometry method provides a rapid and cost-effective tool for cell concentration measurements, as well as cell cycle analysis, which may benefit smaller research laboratories, where flow cytometer is not readily available.

2. INTRODUCTION

Cell cycle analysis is one of the key components of clinical diagnosis and biomedical research. Cell cycle analysis distinguishes cells in different phases of cell cycle and is often used to determine cellular response to drugs and biological stimulations. The assay is based on measuring DNA content of cells to determine if the cells are behaving normally under certain perturbation to the environment such as drug induction. By using a nucleic acid stain, Propidium iodide, one could use fluorescence to accurately analyze cell cycle.

3. CELLOMETER® VISION IMAGING CYTOMETRY

Bright-field (BR) and fluorescent (FL) images

Counting Results

Sample: HCT 116_GFP	Dilution: 1:0
Assay: Annexin V-FITC_Transfection Rate 1	Description: Cell line transfection rate measurement using GFP
Cell: GFP	
Count	240
Mean Size	13.46
Concentration	9.76 x 10 ⁵
Fl Count	56.9%

4. CELL CYCLE DOSE RESPONSE EXPERIMENT

Instrumentation

- Filter Optics Module
- Propidium Iodide Detection
- VB-595-501

Consumable

- 20 µL sample volume
- Concentration ranges: 10⁵ – 10⁷ cells/ml

Drug Induced Cell Cycle Phase Arrest

Experimental Protocol

- Jurkat cells are cultured under 37°C and 5% CO₂
- 4 separate flasks are induced with various concentrations of Aphidicolin, Etoposide, and Nocodazole for 24 h.
- Aphidicolin (1.2, 6, 30 µg/ml)
- Etoposide (0.12, 0.6, 3 µM)
- Nocodazole (.004, .02, .1 µg/ml)

5. CELL CYCLE ANALYSIS AND THEORETICAL RESULTS

Cell Cycle Analysis

- The Jurkat cells induced with the 3 drugs are spun down at 1500 rpm for 5 min
- The supernatant is replaced with 66% ethanol and kept on ice for 15 min
- The cells are again spun down and replaced with 50 µg/ml Propidium iodide
- The cells are allowed to incubate with stain for 30 min before Cellometer Vision and flow cytometry analysis

Aphidicolin

- Aphidicolin induced Jurkat cells should arrest in the G₀/G₁ and S phase

Etoposide

- Etoposide induced Jurkat cells should arrest in the S phase

Nocodazole

- Nocodazole induced Jurkat cells should arrest in the G₂/M phase

6. CELL CYCLE PHASE ARREST RESULTS

Propidium Iodide Fluorescent Image

Cellometer Vision

Flow Cytometer

Aphidicolin 30 µg/ml

More dimmer cells, single copy of DNA (N)


Etoposide 3 µM

Both N and 2N cells, within the sample

Nocodazole 0.1 µg/ml

More brighter cells, duplicate copies of DNA (2N)

8. CELLOMETER VISION AND FLOW CYTOMETER COMPARISON



9. CONCLUSION & ACKNOWLEDGEMENT

We have demonstrated new cell cycle analysis application for the Cellometer imaging cytometry method, which is of great importance to a variety of clinical diagnosis and biomedical research. The ability to rapidly and effectively perform cell population analysis assays, as well as cell concentration and viability measurement, may improve research productivity, especially where flow or laser scanning cytometers are not readily available.

The authors would like to thank the Biotechnology Center at the Merrimack College, North Andover, MA for providing support in this work.