

Cellometer® Auto 1000

Cell Counter

Accurate Cellometer Cell Counting ... Now Even Faster and Easier

- ➔ Improve counting of live/dead cells, individual cells within clusters, and irregular-shaped cells: Proprietary Cellometer software
- ➔ Spend less time counting and more time on research: Concentration and viability in <10 seconds!
- ➔ Convenient stand-alone system with a large, user-friendly touch screen
- ➔ View, print, and save counted cell images, data reports, and cell size histograms: Connect to a network or use a USB drive




Simple, Automated Cell Counting

Improve Speed, Accuracy, and Consistency




Load 20 µl of Sample




Insert Counting Chamber



Select Assay & Click Count



View Image & Results



For more information on the new Cellometer Auto 1000, contact Nexcelom at 978-327-5340 or info@nexcelom.com!



w: www.nexcelom.com
p: 978.327.5340
e: info@nexcelom.com

Cellometer®
Simply Counted

Celigo®
Image Cytometer

Nexcelom products are for RESEARCH USE ONLY and are not approved for diagnostic or therapeutic use.
© Copyright 2014 Nexcelom Bioscience LLC. All Rights Reserved.

Nexcelom
Bioscience

1001197 Rev.C 07/14

